




Instructions for Use




Balance™ Knee OFM1

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	15
ES	Instrucciones para el uso	20
IT	Istruzioni per l'uso	26
NO	Bruksanvisning	31
DA	Brugsanvisning	36
SV	Bruksanvisning	41
EL	Οδηγίες Χρήσης	46
FI	Käyttöohjeet	51
NL	Gebruiksaanwijzing	56
PT	Instruções de Utilização	62
PL	Instrukcja użytkowania	67
TR	Kullanım Talimatları	72
RU	Инструкция по использованию	77
JA	取扱説明書	83
ZH	中文说明书	88
KO	사용 설명서	92
CS	Návod k použití	97
SK	Indikácie použitia	102
RO	Indicații de utilizare	107
HR	Upute za uporabu	112
HU	Használati javallat	117
BG	Указания за употреба	122
SL	Indikácie pre použitie	128
UA	Інструкція з використання	133

DESCRIPTION

The device is a four-bar polycentric mechanical knee. It is designed to facilitate single cadence ambulation on level ground and small slopes.

The device offers a manual locking for improved security which can be activated optionally, aimed for users that require a manual locking feature during rehabilitation.

Extension of the knee is aided by a spring, which is compressed during knee flexion.

The device has a distal tube connector and is supplied with one of these four proximal adapter options (**Fig. 1**). Only replace the adapter with one of these adapter options:

- Pyramid Adapter (A)
- 36 mm Thread (B) for use with 3-Prong Adapter
- Loop Adapter (C)
- IKF Adapter (D)

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low impact use, e.g., gentle walking.

The weight limit for the device is 136kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

36 mm Thread for use with 3-Prong Socket Adapter

Please refer to 3-Prong Socket Adapter Instructions for use.

After lamination, screw the threaded adapter (**Fig. 1: B**) completely into the 3-prong adapter and tighten the lock screw to 10 Nm. After adjustments, the lock screw must be secured with medium strength threadlocker.

Caution: When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the device.

Caution: When carrying out any later adjustments, do not unscrew the threaded adapter more than one full turn. Threaded adapter might break.

Loop Adapter (Fig. 1)

Lamination of adapter to socket:

- Create a dummy disc to replace the adapter disc during lamination. Fix it to the loop adapter (E).
- Laminate the loop adapter (E) to the socket.
- Protect adapter threads during lamination.
- Remove the dummy disc and attach adapter disc (F) before attaching the loop adapter to the knee.

Caution: Fill the central thread of the loop adapter completely with Plastilin, so that no resin can run into the thread during lamination and the central screw can later be completely screwed in.

Caution: Check that the fine serrations on the inner surfaces of the adjustment discs face towards each other.

All Adapters (Fig. 1)

1. Fit the Washers (L and M) onto the Central Screw (I, J or K),
 - match the conical surfaces together.
 - match the fine serrations surfaces together.
2. Mount the adapter on the knee with the Central Screw, using medium strength threadlocker.
3. Tighten the screw to 35 Nm.

Warning: Screws for other adapters may have different lengths. Do not mix the central screws. Check screw for full thread engagement.

Caution: The IKF adapter has reference marks A and P for Anterior and Posterior (G). Make sure to position correctly.

Mounting the Manual Release

The manual release should only be mounted following full prosthesis assembly, length set up, and alignment.

For laminated sockets

1. Laminate the Screw Insert in the desired position.
2. Fix the Cable Guide in the laminate.
3. Mount the Dual Function Lever on the Screw Insert.
4. Secure the Lever Fixation Screw with a medium strength threadlocker.
5. Guide the Nylon Cord through the Cable Guide.
6. Set the Dual Function Lever in the lower position.
7. Snap the Release Lever (**Fig. 7: A**) into the device.
8. Fix the Nylon Cord to the Release Lever on the device.

For plastic sockets

1. Drill a hole in the socket in the desired position and fix the Screw Insert from the inside of the socket.
2. Mount the Dual Function Lever on the Screw Insert.
3. Secure the Lever Fixation Screw with a medium strength threadlocker.
4. Use the Guide Clips to determine where the Cable Guide will run.
5. Fix the Cable Guide in the Guide Clips.
6. Guide the Nylon Cord through the Cable Guide.
7. Set the Dual Function Lever in the lower position.
8. Snap the Release Lever (**Fig. 7: A**) into the device.
9. Fix the Nylon Cord to the Release Lever on the device.

Caution: If the device is being used with the locking function, neither the Release Lever (**Fig. 7: A**) nor the Nylon Cord may be impeded by the cosmesis and they must be able to move freely.

Warning: If readjustment work is required, secure the Nylon Cord so that it cannot catch in grinding tools.

Note: If the patient only wants to use the Dual Function Lever for releasing the knee joint in certain situations (e. g., sitting down), the enclosed M4 Allen Screws can be screwed into the slots in the Dual Function Lever. This stops the lever from slipping up beyond the taper and thus permanently inactivating the Release Lever (**Fig. 7: A**).

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 2)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through midtube (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through midtube (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 3**).

Tighten the screws with the following torque:

- Central Screw (**Fig. 1: I, J, or K**): 35 Nm
- Tube Clamp Screw (**Fig. 4: A**): 16 Nm

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

Depending on which device setting is used, ensure that the patient is familiar with the operation of the manual lock release function, or the geometric locking function.

Note: The device is configured as a conventional locking knee at the factory and is delivered in this setting. To use the device as a freely movable joint, please refer to the following section:

Converting the Joint From Locking to Mechanical.

Adjustment of stance phase stability (Fig. 5)

Stance phase stability can be adjusted by shifting the knee according to the following:

- To increase stance phase stability; shift knee backwards compared to socket (A).
- To allow for a more dynamic gait: shift knee forward compared to socket (B).

For IKF adapter only: Stance flexion adjustment

The back bumper (**Fig. 1: H**) can be changed to increase or decrease stance flexion:

Back Bumper (with holes)		
Yellow	Soft	More stance phase damping
Red	Medium	Standard
Brown	Hard	Less stance phase damping

To replace a bumper, unscrew the central screw. Replace the bumper (holes of back bumper facing up) and reassemble according to assembly instructions above (torque to 35 Nm).

Swing control

Adjusting the Extension Assist Spring (**Fig. 6**):

Turn the Extension Spring Adjustment Screw (A) to the:

- right to increase the action of the Extension Spring.
- left to decrease the action of the Extension Spring.

Adjusting axle friction (Fig. 7):

Turn the Friction Screw (C) to the:

- right to increase friction on the axle.
- left to decrease friction on the axle.

Note: Swing control is working only when using the device as a freely movable joint with integrated Extension Spring (refer to the following section: **Converting the Joint From Locking to Mechanical**).

Converting the Joint (Fig. 7)

From Locking to Mechanical

- Push the Release Lever (A) upwards and hold it in this position.
- Turn the two Release Screws (B) adjacent to the Release Lever (A) finger-tight to the right to secure it in this position.
- Tighten the Release Screws (B) to 5 Nm and always secure them with a medium strength threadlocker when finishing the prosthesis.

Caution: Failure in above mentioned conversion when tightening the Release Screws (B) results in loss of function of the device and/or makes it impossible to release the joint.

Caution: If the device is being used with the freely moving function, the Release Screws (B) must be tightened to 5 Nm and secured with a medium strength threadlocker.

From Mechanical to Locking

- Conversion to a fixed knee is performed in the reverse order by turning the two Release Screws (B) to the left until the Release Lever (A) moves freely.
- Secure the Release Screws (B) again with a medium strength threadlocker and check that the Release Lever (A) is freely movable.

Caution: After every conversion, carry out a functional check before the dynamic alignment.

Caution: If the device is being used with the locking function, the Release Screws (B) must not impede the function of the Release Lever (A) and must be secured with a medium strength threadlocker to rule out any inadvertent loosening and incorrect function.

Note: If the device is being used with the locking function, a certain amount of play may develop. This play can be eliminated with the help of the Extension Stop Adjustment Screw (D) which is

turned right. This moves the Extension Stop, which compensates for the play. If the Release Lever (A) does not click into place easily, correct the position of the Extension Stop with the help of the Extension Stop Adjustment Screw (D), which you turn to the left.

Caution: When readjusting the Extension Stop with the Extension Stop Adjustment Screw (D), make sure that the device functions or locks properly and carry out a functional check before the dynamic alignment.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -10°C to 40°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P6 - 136 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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