
Instructions for Use


---

BALANCE™ KNEE OFM2


|    |                             |     |
|----|-----------------------------|-----|
|    |                             | 3   |
| EN | Instructions for Use        | 4   |
| DE | Gebrauchsanweisung          | 9   |
| FR | Notice d'utilisation        | 15  |
| ES | Instrucciones para el uso   | 21  |
| IT | Istruzioni per l'uso        | 27  |
| NO | Bruksanvisning              | 33  |
| DA | Brugsanvisning              | 38  |
| SV | Bruksanvisning              | 44  |
| EL | Οδηγίες Χρήσης              | 50  |
| FI | Käyttöohjeet                | 56  |
| NL | Gebruiksaanwijzing          | 61  |
| PT | Instruções de Utilização    | 67  |
| PL | Instrukcja użytkowania      | 73  |
| TR | Kullanım Talimatları        | 79  |
| RU | Инструкция по использованию | 84  |
| JA | 取扱説明書                       | 90  |
| ZH | 中文说明书                       | 95  |
| KO | 사용 설명서                      | 100 |


1


2


3

4

5

### DESCRIPTION

The device is a weight-activated friction brake knee intended to provide stability in stance through alignment or via a weight activated friction brake.

The device has an internal extension assist, a spring-loaded mechanism for controlling swing properties and returning knee to full extension at end of swing phase. This extension assist is adjustable to individual patients.

The device has a Distal Tube Connector and a Proximal Integrated Male Pyramid Adapter.

### INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low impact use, e.g., gentle walking.

The weight limit for the device is 125 kg.

### GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

**Warning:** Avoid placing hands or fingers near moving joints.

**Caution:** Do not adjust other screws than described in these instructions.

The device is for single patient use.

### ALIGNMENT INSTRUCTIONS

#### ***Bench Alignment (Fig. 1)***

#### **Alignment Goal**

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through midtube (A)
- fall at the 1/3 mark on the inside of the foot cover.

**Note:** Prioritize knee alignment over foot alignment if there is a mismatch.

## Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through midtube (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

**Warning:** After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

**Caution:** Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 2**).

Tighten the screws with the following torque:

- Tube Clamp Screw (**Fig. 3**): 16 Nm

## Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

## Dynamic Alignment

Ensure the patient is familiar with the functioning of the device. Depending on which device setting is used, ensure that the patient is familiar with the operation of the manual lock release function, or the weight-activated brake function.

**Note:** The device is configured as a conventional locking knee at the factory and is delivered in this setting. To use the device as a freely movable joint, please refer to the following section: **Converting the Joint - From Locking to Mechanical**.

## Stance phase stability (Fig. 4):

Start carefully with the device's brake function with the patient at the dynamic alignment and check whether this basic setting matches the patient's needs. If not, please proceed as follows:

## Adjustment using the Load Regulating Screw (4A)

The Load Regulating Screw can be used to adjust:

- the point at which the brake acts when (weight-) loading starts,
- and the response of the brake during the dynamic alignment.

Variant 1 - Brake effect shall be induced later:

Turn the Load Regulating Screw (4A) to the right which increases the weight needed by the patient to initiate the brake function.

Variant 2 - Brake effect shall be induced earlier:

Turn the Load Regulating Screw (4A) to the left which reduces the weight needed by the patient to initiate the brake function.

**Note:** Always adjust the Load Regulating Screw (4A) carefully. Note that you can adjust it against the slight resistance of the threadlocker.

Readjusting the break unit after prolonged use

After a prolonged period of use, play may develop in the device which can be readjusted using the Service Screw (4B).

Turn the Service Screw by small steps to the left.

Carefully check whether the brake function of the device now suits the patient's individual needs. If not, readjust the Service Screw again until the correct setting has been found.

**Note:** Always adjust the Service Screw (4B) carefully. Note that you can adjust it against the slight resistance of the threadlocker.

**Caution:** Flexion and extension must be possible at all settings.

Restoring the basic setting

If necessary, the basic setting of the device can be reset as follows:

Step 1: Carefully turn the Load Regulating Screw (4A) to the right until a slight resistance is felt. Then turn the Load Regulating Screw two turns back to the left.

Step 2: Carefully turn the Service Screw (4B) to the right until a slight resistance is felt. Then turn the Service Screw three turns back to the left. It must now be possible to flex the joint again.

**Caution:** Flexion and extension must be possible at all settings.

## Swing control

Adjusting the Extension Assist Spring (**Fig. 5**):

Turn the Extension Spring Adjustment Screw (A) to the:

- right to increase the action of the Extension Spring.
- left to decrease the action of the Extension Spring.

**Note:** Swing control is working only when using the device as a freely movable joint with integrated Extension Spring (refer to the following section: **Converting the Joint - From Locking to Mechanical**).

## ASSEMBLY INSTRUCTIONS

For laminated sockets

1. Laminate the Screw Insert in the desired position.
2. Fix the Cable Guide in the laminate.
3. Mount the Dual Function Lever on the Screw Insert.
4. Secure the Lever Fixation Screw with a medium strength threadlocker.
5. Guide the Nylon Cord through the Cable Guide.
6. Set the Dual Function Lever in the lower position.
7. Snap the Release Lever (**Fig. 4: C**) into the device.
8. Fix the Nylon Cord to the Release Lever on the device.

For plastic sockets

1. Drill a hole in the socket in the desired position and fix the Screw Insert from the inside of the socket.
2. Mount the Dual Function Lever on the Screw Insert.
3. Secure the Lever Fixation Screw with a medium strength threadlocker.
4. Use the Guide Clips to determine where the Cable Guide will run.

5. Fix the Cable Guide in the Guide Clips.
6. Guide the Nylon Cord through the Cable Guide.
7. Set the Dual Function Lever in the lower position.
8. Snap the Release Lever (**Fig. 4: C**) into the device.
9. Fix the Nylon Cord to the Release Lever on the device.

**Caution:** If the device is being used with the locking function, neither the Release Lever (**Fig. 4: C**) nor the Nylon Cord may be impeded by the cosmesis and they must be able to move freely.

**Warning:** If readjustment work is required, secure the Nylon Cord so that it cannot catch in grinding tools.

**Note:** If the patient only wants to use the Dual Function Lever for releasing the knee joint in certain situations (e. g., sitting down), the enclosed M4 Allen Screws can be screwed into the slots in the Dual Function Lever. This stops the lever from slipping up beyond the taper and thus permanently inactivating the Release Lever (**Fig. 4: C**).

### ***Converting the Joint (Fig. 4)***

#### **From Locking to Mechanical**

Push the Release Lever (C) upwards and hold it in this position.

Turn the Release Screw (D) adjacent to the Release Lever (C) finger-tight to the right to secure it in this position.

Tighten the Release Screw (D) to 5 Nm and always secure it with a medium strength threadlocker when finishing the prosthesis.

**Caution:** Failure in above mentioned conversion when tightening the Release Screw (D) results in loss of function of the device and/or makes it impossible to release the joint.

**Caution:** If the device is being used with the freely moving function, the Release Screw (D) must be tightened to 5 Nm and secured with a medium strength threadlocker.

#### **From Mechanical to Locking**

- Conversion to a fixed knee is performed in the reverse order by turning the Release Screw (D) to the left until the Release Lever (C) moves freely.
- Secure the Release Screw (D) again with a medium strength threadlocker and check that the Release Lever (C) is freely movable.

**Caution:** After every conversion, carry out a functional check before the dynamic alignment.

**Caution:** If the device is being used with the locking function, the Release Screw (D) must not impede the function of the Release Lever (C) and must be secured with a medium strength threadlocker to rule out any inadvertent loosening and incorrect function.

**Note:** If the device is being used with the locking function, a certain amount of play may develop. This play can be eliminated with the help of the Extension Stop Adjustment Screw (E) which is turned right. This moves the Extension Stop, which compensates for the play.

If the Release Lever (C) does not click into place easily, correct the position of the Extension Stop with the help of the Extension Stop Adjustment Screw (E), which you turn to the left.

**Caution:** When readjusting the Extension Stop with the Extension Stop Adjustment Screw (E), make sure that the device functions or locks properly and carry out a functional check before the dynamic alignment.

## USAGE

### *Cleaning and care*

Wipe the device with a soft cloth. Do not use solvents.

### *Environmental Conditions*

The device should not come in contact with fresh water, salt water or chlorinated water.

**Caution:** Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15 °C to 50 °C.

## MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

**Caution:** Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

**Caution:** Do not use talcum powder for lubrication.

## REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

### *Compliance*

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P6 - 125 kg \*)

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use  
see manufacturer's written instructions on  
intended use!

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