



## Instructions for Use













# Cheetah® Xplore Junior

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	16
ES	Instrucciones para el uso	21
IT	Istruzioni per l'uso	26
NO	Bruksanvisning	30
DA	Brugsanvisning	35
SV	Bruksanvisning	40
EL	Οδηγίες Χρήσης	44
FI	Käyttöohjeet	50
NL	Gebruiksaanwijzing	54
PT	Instruções de Utilização	59
PL	Instrukcja użytkowania	64
TR	Kullanım Talimatları	69
RU	Инструкция по использованию	74
JA	取扱説明書	79
ZH	中文说明书	84
KO	사용 설명서	88
CS	Návod k použití	92
SK	Indikácie použitia	97
RO	Indicații de utilizare	102
HR	Upute za uporabu	107
HU	Használati javallat	111
BG	Указания за употреба	116
SL	Indikacije pre poročilo	122
UA	Інструкція з використання	126

## DESCRIPTION

The device is a prosthetic foot attached to the posterior wall of the socket. It is composed of a Cheetah Junior inspired toe blade and a sole blade secured using heel bolts.

The device consists of the following components:

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This device must be used with a Foot Cover and a Spectra Sock.

## INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is intended for children and adolescent patients.

## GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** The device may break or delaminate if sudden high impact side bending and/or twisting occurs.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

**Caution:** Epoxy can cause allergic reactions when in direct contact with body.

## REQUIRED COMPONENTS

The following components are required to use the device:

- Lamination plate
- Attachment hardware

## DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

**Warning:** Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	15–24	25–34	35–45	46–55
Weight lbs	35–53	54–75	76–99	100–121
Category	1	2	3	4

## ASSEMBLY INSTRUCTIONS

### **Foot Cover and Spectra Sock (Fig. 2)**

**Caution:** To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

**Note:** The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

### **Prosthesis**

Assemble prosthesis with applicable devices.

**Warning:** Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

**Warning:** Ensure proper attachment by following the applicable device assembly instructions.

### **Sole Plate Attachment**

- Grind the sole plate to appropriate size, using the grinding template included in the foot kit (**Fig. 4**). The size of the foot is 2 cm shorter when used with the Foot Cover. This needs to be taken into account when grinding the sole plate.
- Use bolts (included in foot kit) to attach sole plate to foot module. Apply Loctite 410 to all threads and torque to following specifications: 14 Nm
- Follow alignment below (**Fig. 3**).

**Note:** If foot will be used without shoes, proper sole material needs to be bonded to the sole plate for protection of the foot module.

### **Lamination**

The foot may be connected with either Lamination Connector or direct lamination.

#### **Lamination Connector**

The Lamination Connector must be laminated directly into the posterior wall of the socket.

**Warning:** Before starting make sure to orient the Lamination Connector the correct way. The shoulders of the inserts (marked with arrow) should be oriented towards the socket (**Fig. 5**).

1. Put the Lamination Connector in the middle third of the socket to make height adjustment easier.
2. Make sure the socket surface is flat. A flat surface below the Lamination Connector distributes the load to the connected foot.
3. Use a friction pad between the Lamination Connector and the foot.
4. Apply Loctite 410 and torque bolts to 28 Nm for bolt size 5/16".

#### **Cutting Foot Module**

- Check for clearance in swing phase.
- Trim foot module to the top end of the laminating connector.
- Cut foot module using a hack saw or band saw.

## **Direct Lamination**

Before you start you need a well fitting conventional prosthesis. Either an existing socket or alternatively a prosthesis with a check socket on conventional prosthesis can be used.

### **1. Laminate the Socket (Fig. 6)**

- a. Laminate the socket using carbon braid. Add extra layers of carbon in the upper posterior surface where the foot will be attached.
- b. Roughen the upper posterior surface to improve the adhesion to the socket when the foot module is attached (shown by grey area).

### **2. Trace around the Original Walking Prosthesis (Fig. 7)**

- a. Ensure that the prosthetic foot is resting firmly on the shoes insole and fasten the laces tightly.
- b. Lay the original walking prosthesis over the grid on the A1 alignment sheet (included)
- c. Mark the lowest points of the anterior and posterior brims as a reference.
- d. With the shoe on, trace a line around the entire prosthesis. Then remove the shoe and trace the outline of the foot cover, ensuring the prosthesis remains aligned with the original tracing.

### **3. Fill in the Patient Information Table on the Alignment Sheet (Fig. 8)**

- a. The socket adduction can be observed by looking at the original walking prosthesis in the frontal plane.
- b. Use the angle guide on the A1 alignment sheet to help measure the foot rotation in relation to the socket.

### **4. Cut the Pylon to Length (Fig. 9)**

- a. Align the new socket and the foot with the tracing of the original prosthesis.
- b. Mark the pylon at the height of the posterior brim of the socket using the drawing of the original prosthesis.
- c. Cut the pylon at this mark and round off the end, removing any sharp edges.

### **5. Mark the Length of the Attachment Surface on the Foot Pylon (Fig. 10)**

- a. Mark the pylon at 90 mm from the top (cut end).
- b. Roughen all the surfaces of the pylon above this mark (shown by grey area).

### **6. Align the Foot and the Socket with the Tracing of the Original Walking Prosthesis (Fig. 11)**

- a. The sole should be lined up with the outline of the foot cover sole.
- b. The new socket should be positioned over the tracing of the original socket with additional 5-7° of knee extension in the socket.
- c. Use a foam spacer/s (or other spacer) placed between the pylon and the socket to maintain the correct extension. Spacers should be placed at the end of the roughened attachment surface. Ensure the minimum clearance of 5 mm is maintained between the distal posterior area of the socket and the foot pylon.
- d. Any spacer used at the top of the pylon can have maximum thickness of 5 mm.
- e. Fasten the foam to the foot using tape.
- f. Mark a clear height on the socket at height of a foam spacer on the pylon.

### **7. Align the Socket according to the Patient Information (Fig. 12)**

- a. Put the shoe on the foot and place it upright in the work surface.
- b. Hold the socket in position on the pylon according to:
  - The height
  - The recorded socket adduction
  - The recorded foot rotation in relation to socket (using angle guides on A1 sheet)
- c. The foam spacers will hold the socket in the correct extension.

### **8. Mark the Attachment Surface onto the Socket (Fig. 13)**

- a. Holding the socket in the desired position, trace around the end of the pylon onto the socket above the height mark.

### 9. Fix the Foot to the Socket (**Fig. 14**)

- a. Use Composite 1 adhesive that is provided with the foot kit, read all information on the adhesive packaging.  
**Note:** Adhesive working time is 1 minute.
- b. Cut the tip of the adhesive mixer at 16 mm from the end.
- c. Apply the adhesive generously onto the marked attachment surface on the socket.
- d. Attach the foot onto the marked attachment area and make any small adjustment to ensure the alignment is correct.
- e. Smooth the surface of any adhesive that has been pushed out around the outside of the pylon.
- f. Hold the socket in position for 1 minute, then leave for a further 5 minutes for the adhesive to set.

### 10. Apply Adhesive around the outside of the Pylon (**Fig. 15**)

- a. Apply adhesive generously onto the outside roughed surface of the pylon and spread outwards onto the roughened surface of the socket.
- b. Smooth the finish of the adhesive and allow to set for 5 minutes.

### 11. Test Alignment before Final Lamination (**Fig. 16**)

- a. Remove any foam spacers and tape before testing the alignment on the patient.
- b. If needed, minor adjustments can be made by shaping the underside of the sole.
- c. If the alignment needs further adjustment, grind/cut the foot off the socket and make the desired changes before repeating steps 10 and 11.
- d. If the alignment is good, proceed to laminate over the attachment area onto the socket using carbon braid and a hand lay up.

## ALIGNMENT INSTRUCTIONS

### ***Bench Alignment (Fig. 3)***

#### **Sagittal plane**

- Determine vertical socket position in alignment jig.
- Adjust to appropriate heel height.
- Introduce appropriate socket angles flexion/extension and abduction/adduction.
- The load line should fall at the junction of the posterior and the middle third portions as shown in **Fig. 3**.

#### **Coronal Plane**

Introduce appropriate socket abduction/adduction using alignment jig.  
Determine appropriate internal or external socket rotation.

## USAGE

### ***Cleaning and Care***

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

### ***Environmental Conditions***

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

Dry with a cloth after contact with fresh water or humidity.

## MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

### **Noise from Foot**

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

### **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### **LIABILITY**

Össur devices are designed and verified to be safe and compatible in combination with each other, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

### **Compliance**

This device has been tested according to ISO 10328 standard to two million load cycles. Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - Label		
Category	Weight (Kg) Extreme Impact	Label Text
1	24	ISO 10328-P3*-24 kg
2	34	ISO 10328-P3*-34 kg
3	45	ISO 10328-P3-45 kg
4	55	ISO 10328-P3-55 kg

\*Actual test loads reflect maximum body mass

ISO 10328 - "P" - "m"kg \*)

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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