



## Instructions for Use

# Flex-Foot Assure®

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	16
IT	Istruzioni per l'uso	21
NO	Bruksanvisning	25
DA	Brugsanvisning	29
SV	Bruksanvisning	33
EL	Οδηγες χρήσης	38
FI	Käyttöohjeet	42
NL	Gebruiksaanwijzing	46
PT	Instruções de Utilização	50
PL	Instrukcja użytkowania	55
TR	Kullanım Talimatları	59
RU	Инструкция по использованию	63
JA	取扱説明書	67
ZH	中文说明书	72
KO	사용 설명서	75
CS	Návod k použití	79
SK	Indikácie použitia	83
RO	Indicații de utilizare	87
HR	Upute za uporabu	92
HU	Használati javallat	96
BG	Указания за употреба	100
SL	Indikacije pre poročilo	104
UA	Інструкція з використання	109

## DESCRIPTION

The device is a prosthetic foot with carbon heel and keel. Designed specifically for slower walkers, Flex-Foot Assure incorporates an active heel combined with a full-length keel that work together to protect the vulnerable sound limb.

The device consists of the following components (**Fig. 1**):

- A. Foot Module
- B. Foam Wedge
- C. 30 mm Fiberglass Tube Kit
- D. Tube Insert for SF5-SF6
- E. 30 mm Alignment Fitting Tool
- F. Male Pyramid
- G. Attachment Bolts

This device must be used with a Foot Cover and a Spectra Sock.

## INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low impact use, e.g., gentle walking.

## GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

## DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

**Warning:** Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	45-61	62-80	81-97	98-114	115-136
Weight lbs	99-135	136-176	177-214	215-251	252-300
Size 19-20	1	2	N/A	N/A	N/A
Size 21-22	1	2	N/A	N/A	N/A
Size 23-24	1	2	3	4	N/A
Size 25-26	1	2	3	4	5

Weight kg	45-61	62-80	81-97	98-114	115-136
Weight lbs	99-135	136-176	177-214	215-251	252-300
Size 27-28	1	2	3	4	5
Size 29-30	N/A	2	3	4	5

## ASSEMBLY INSTRUCTIONS

### *Foot Cover and Spectra Sock (Fig. 2)*

**Caution:** To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

**Note:** The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

### **If required doff the Foot Cover as follows:**

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

### **Bolt Assembly**

To attach a Pyramid or a Tube, thread each bolt through washer and Foot Module into Pyramid/Tube. Apply Loctite 410 to the bolt and torque according to the relevant instructions. When using a Tube, place friction pad between foot and Tube.

### **Prosthesis**

Assemble prosthesis with applicable devices.

**Warning:** Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Remove the protective film on the pyramid after fitting.

## ALIGNMENT INSTRUCTIONS

### **Heel Wedge**

- Provides smooth roll over and improved cushioning at heel strike.
- The heel wedge is an essential component for Flex-Foot Assure function.

**Note:** the heel wedge should always be used with the foot.

### **Bench Alignment (Fig. 3)**

#### **Alignment Goal**

Alignment reference line (B) should:

- pass through midpoint of socket at either the patellar tendon level or at the ischial tuberosity level (D).
- fall at the 1/3 mark on the inside of the foot cover.

**Note:** Prioritize knee alignment over foot alignment if there is a mismatch.

#### **Alignment Instructions**

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect either the socket or the knee to the foot and establish the correct knee center height.

3. If using a prosthetic knee: Position knee according to knee alignment instructions
4. On the lateral side of the socket, make a first mark (alignment reference point) at the midpoint of the socket at either patellar tendon level or the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through the alignment reference point (D).
6. Set the correct socket angles for flexion/extension and abduction/adduction.
7. If using a prosthetic knee: Use the applicable adapters to connect the socket to the prosthesis.

### **Static Alignment**

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

### **Dynamic Alignment**

Ensure the patient is familiar with the functioning of the device.

The heel to toe action can be influenced by:

- Heel Stiffness
- Anterior-Posterior positioning of device
- Dorsi-Plantarflexion
- Shoe Characteristics

### **Consider the following actions if needed:**

Symptoms

- Device comes to flat position too early (patient feels he/she is sinking into a hole)
- Climbing over the toe requires extra energy
- Toe feels too stiff
- Knee hyperextends

Action

- Check Wedge instructions (if applicable)
- Shift socket anterior (or device posterior)
- Consider dorsiflexion
- Check heel of the shoe and shoe performance

Symptoms

- Rapid heel to toe movement
- Poor control over prosthesis at initial contact
- Minimal energy return feeling
- Too little push off from the toe
- Knee becomes unstable

Action

- Check Wedge instructions (if applicable)
- Shift socket posterior (or device anterior)
- Consider plantarflexion
- Check heel of the shoe and shoe performance

## **USAGE**

### **Cleaning and care**

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

### **Environmental Conditions**

The device is Waterproof if used with Tube.

The device is Weatherproof if used with Pyramid.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps. It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e.g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

## **MAINTENANCE**

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

### **Noise from Foot**

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

## **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:


- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

## **Compliance**

This device has been tested according to ISO 10328 standard to two million load cycles. Depending on patient activity this may correspond to 2-3 years of use.

<b>ISO 10328 - Label</b>		
<b>Category</b>	<b>Weight (Kg) Low Impact</b>	<b>Label Text</b>
0	44	ISO 10328-P3*-44 kg
1	61	ISO 10328-P3*-61 kg
2	80	ISO 10328-P3-80 kg
3	97	ISO 10328-P5-97 kg

ISO 10328 - Label		
Category	Weight (Kg) Low Impact	Label Text
4	114	ISO 10328-P6-114 kg
5	136	ISO 10328-P6-136 kg
*Actual test loads reflect maximum body mass		

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