
Instructions for Use

FLEX-FOOT[®] BALANCE






Product Number:

FBP0xyz(BR)

FBPUxyz(BR)

| | | |
|----|-----------------------------|----|
| | | 3 |
| EN | Instructions for Use | 4 |
| DE | Gebrauchsanweisung | 6 |
| FR | Notice d'utilisation | 9 |
| ES | Instrucciones para el uso | 11 |
| IT | Istruzioni per l'uso | 14 |
| DA | Brugsanvisning | 16 |
| SV | Bruksanvisning | 18 |
| EL | Οδηγίες χρήσης | 20 |
| NL | Gebruiksaanwijzing | 23 |
| PT | Instruções de Utilização | 26 |
| PL | Instrukcja użytkowania | 29 |
| CZ | Návod k použití | 32 |
| TR | Kullanım Talimatları | 35 |
| RU | Инструкция по использованию | 37 |
| JA | 取扱説明書 | 40 |
| ZH | 中文说明书 | 42 |
| KO | 사용 설명서 | 44 |

ENGLISH

FLEX-FOOT® BALANCE FOOT SYSTEM (Figure 1)

- Build height: 56-62 mm (2 3/16" - 2 7/16")
- Maximum user weight: 136 kg (299 lbs)
- Available in 10 mm (3/8") heel height
- Available in sizes 21-28 cm

NOTE: Select foot category according to impact level.

Selection chart is provided in the catalog and on website www.ossur.com

BENCH ALIGNMENT (Figure 2 & 3)

- Fit foot with selected cover.
- Adjust to appropriate heel height (using the shoe).
- Introduce appropriate socket angles flexion/ extension and abduction/ adduction.
- Divide the foot cover length into 3 equal portions.
- The load line should fall at the junction of the posterior and the middle third portions as shown in the diagram.

DYNAMIC ALIGNMENT

The heel stores energy after initial contact, slowly releasing it at midstance. The heel to toe action can be influenced by

- A-P positioning of foot
- Dorsi-Plantarflexion

CATEGORY SELECTION CHART

Please refer to the selection charts below to determine the appropriate stiffness required according to Össur recommendations.

| WEIGHT KG | 36-44 | 45-61 | 62-80 | 81-97 | 98-114 | 115-136 |
|------------|-------|--------|---------|---------|---------|---------|
| WEIGHT LBS | 79-97 | 98-134 | 135-176 | 177-214 | 215-251 | 252-300 |
| 21 | 0 | 1 | N/A | N/A | N/A | N/A |
| 22 | 0 | 1 | N/A | N/A | N/A | N/A |
| 23 | 0 | 1 | 2 | N/A | N/A | N/A |
| 24 | 0 | 1 | 2 | N/A | N/A | N/A |
| 25 | N/A | 1 | 2 | 3 | 4 | 5 |
| 26 | N/A | 1 | 2 | 3 | 4 | 5 |
| 27 | N/A | N/A | 2 | 3 | 4 | 5 |
| 28 | N/A | N/A | 2 | 3 | 4 | 5 |

HEEL STIFFNESS

Heel too soft

Symptoms

- Foot will come to flat position too early (amputee feels he/she is sinking into a hole).
- Extra energy is required to climb up over the toe.
- Toe will feel too stiff.
- Knee may hyperextend.

Actions

- Shift socket anteriorly.
- Dorsi-flex the foot.

Heel too hard

Symptoms

- Rapid heel to toe movement.
- Amputee has poor control of his prosthesis.
- Toe will feel too stiff - minimal energy return.
- Knee may become unstable.

Actions

- Shift socket posteriorly.
- Plantar-flex the foot.

PROTECTIVE FILM ON PYRAMID

The protective film protects the pyramid from scratching. Remove the film when the fitting is complete.

SHOE HORN


When removing and installing the cover use shoehorn to avoid damage to foot or cover (**Figure 4 & 5**).

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

COMPLIANCE

This component has been tested according to ISO 10328 standard to two million load cycles. Depending on the amputee's activity this corresponds to a duration of use of two to three years. We recommend carrying out regular yearly safety checks

In the standard mentioned, test levels (P) are assigned to a certain maximal body masses (m in kg). In some cases, which are marked with, no test level is assigned to the product related maximal body mass. In these cases, the test loads have been

adapted adequately on the basis of the specified load level.

| ISO 10328 - label | | |
|-------------------|-------------|---------------------|
| Category | Weight (Kg) | Label text |
| 0 | 44 | ISO 10328-P3-44-kg |
| 1 | 61 | ISO 10328-P3-61-kg |
| 2 | 80 | ISO 10328-P4-80-kg |
| 3 | 97 | ISO 10328-P5-97-kg |
| 4 | 114 | ISO 10328-P5-114-kg |
| 5 | 136 | ISO 10328-P6-136-kg |

- This manual is intended for use by a certified prosthetist.

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EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If un-usual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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