
Instructions for Use


ÖSSUR®
LIFE WITHOUT LIMITATIONS

Flex-Symes™








		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	16
ES	Instrucciones para el uso	21
IT	Istruzioni per l'uso	27
NO	Bruksanvisning	32
DA	Brugsanvisning	38
SV	Bruksanvisning	43
EL	Οδηγίες Χρήσης	48
FI	Käyttöohjeet	54
NL	Gebruiksaanwijzing	59
PT	Instruções de Utilização	65
PL	Instrukcja użytkowania	71
TR	Kullanım Talimatları	76
RU	Инструкция по использованию	81
JA	取扱説明書	87
ZH	中文说明书	93
KO	사용 설명서	98
CS	Návod k použití	103
SK	Indikácie použitia	108
RO	Indicații de utilizare	114
HR	Upute za uporabu	120
HU	Használati javallat	125
BG	Указания за употреба	131
SL	Indikácie pre použitie	137
UA	Інструкція з використання	142

1

2

DESCRIPTION

The device is a prosthetic foot. The Flex-Symes foot design is a carbon foot module connected to an actively deflecting carbon heel with alignment adapter.

The device consists of the following components (**Fig. 1**):

- A. Foot Module
- B. Heel
- C. Heel Fastener Hardware Kit
- D. Socket Adapter
- E. Alignment Disc Kit
- F. Central Bolt
- G. Anti-Rotation Bolt Kit
- H. Heel Wedge Kit

This device must be used with a Foot Cover and a Spectra Sock.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

Warning: Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-356
Low Impact Level										
Size 22-24	1	1	2	3	4	5	6	N/A	N/A	N/A
Size 25-27	1	1	2	3	4	5	6	7	8	9
Size 28-30	3	3	3	3	4	5	6	7	8	9

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	221-256	257-287	288-324	325-356
Moderate Impact Level										
Size 22-24	1	2	3	4	5	6	N/A	N/A	N/A	N/A
Size 25-27	1	2	3	4	5	6	7	8	9	N/A
Size 28-30	3	3	3	4	5	6	7	8	9	N/A
High Impact Level										
Size 22-24	2	3	4	5	6	N/A	N/A	N/A	N/A	N/A
Size 25-27	2	3	4	5	6	7	8	9	N/A	N/A
Size 28-30	3	3	4	5	6	7	8	9	N/A	N/A

ASSEMBLY INSTRUCTIONS

Lamination Technique

- Mark alignment position after dynamic alignment and remove Foot Module from socket.
Note: If Alignment Wedges are straight, both wedges may be removed to decrease clearance and complexity! This will lower the height by 6 mm.
- Invert the socket and clamp securely in the lamination jig.
- Clean the groove of the aluminum socket adapter with acetone or similar to ensure maximum bonding of the fibers to the aluminum.
- Roughen attachment surfaces and seal center hole (by using plasticine or similar).
- Apply appropriate carbon lay-up (see. Table below). Tie off the overhanging braids around the base groove using a heavy thread. Ensure that the tie stays tight.
- Finish with a stockinette for cosmetic finish, if desired.
- Pull a PVA bag and laminate with acrylic or epoxy resin.
- Tape with PVC-tape around the distal end.
- Tie off excess resin before it reaches the gel state. Allow resin to harden.
- Grind off the excess hardened resin around the edge of the aluminum socket adapter. Expose the attachment surface. Care must be taken not to grind the aluminum piece itself!

Sizes	Category	Layers of Carbon Fiber
22 - 24	1 – 6	Minimum 4*
25 - 30	1 – 6	Minimum 6*
25 - 30	7 – 9	Minimum 8*

*The appropriate number of layers depends on user activity and weight and must be evaluated by the CPO. It is recommended to use heavy carbon fiber braids.

Assembly (Fig. 4, 5 & 6)

- Position circular Alignment Wedges according to Dynamic Alignment.
- Assemble and torque Central Bolt using Loctite 410 (**Fig. 4**).
 - Category 1-6: Torque to 40 Nm
 - Category 7-9: Torque to 60 N

Note: Applies to all sizes.

- Guided by the positioning hole, drill a hole through the circular Alignment Wedge/s and the socket adapter. Tap threads for the Anti-Rotation Bolt (**Fig. 5**).
 - Sizes 22-27: Drill $\varnothing$ 5.0 mm, Tap M6 x 1.0
 - Sizes 28-30: Drill $\varnothing$ 6.8 mm, Tap M8 x 1.25
- Assemble, using Loctite 410, Anti-Rotation Bolt and Anti-Rotation Washer to secure the alignment. Torque to 10 Nm.
- Position the heel divider in the posterior slot of the splitted heel part (**Fig. 6**). Bolt the heel to the Foot Module and torque to specifications using Loctite 410. Make sure there is no interference with the Central Bolt head when the heels are compressed. To secure the position of the heel divider a drop of instant adhesive can be used on one side.
 - Category 1-6: Bolt diameter 1/4 inches, Torque to 14 Nm
 - Category 7-9: Bolt diameter 5/16 inches, Torque to 28 Nm

Note: Applies to all sizes.

Foot Cover and Spectra Sock (Fig. 2)

Caution: To avoid pinching fingers, always use a shoehorn.

The Foot Cover may need some trimming if the socket does not fit into the foot cover. A foot cover attachment cannot be used for Flex-Symes.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

Note: The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

Prosthesis

Assemble prosthesis with applicable devices.

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Warning: Ensure proper attachment by following the applicable device assembly instructions.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 3)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at either the patellar tendon level or at the ischial tuberosity level (D).
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect either the socket or the knee to the foot and establish the correct knee center height.
3. If using a prosthetic knee: Position knee according to knee alignment instructions
4. On the lateral side of the socket, make a first mark (alignment reference point) at the midpoint of the socket at either patellar tendon level or the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through the alignment reference point (D).
6. Set the correct socket angles for flexion/extension and abduction/adduction.
7. If using a prosthetic knee: Use the applicable adapters to connect the socket to the prosthesis.

Fabrication Procedure (Fig. 7)

- Leave foot assembled to the socket adapter and protect foot with appropriate film/sheet barrier.

- Roughen attachment surface, position socket in the desired alignment and fill the space beneath the socket with epoxy putty.
- Let it harden (3-5 minutes).
- Wrap synthetic casting tape and proceed with dynamic alignment for temporary use only.

Note: Do not use polyurethane foam build-up since this compresses and can lead to bolt or lamination failure.

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Alignment Wedges

Use the alignment options of the circular Alignment Wedges to fine tune the dynamic alignment. The Alignment Wedges provide angular adjustments up to 6 degrees.

Instructions:

- Two 3 mm steel pins are provided in the kit.
- Insert pins into alignment wedges.
- Twist Alignment Wedges against each other to fine tune alignment.

Ensure the patient is familiar with the functioning of the device.

The heel to toe action can be influenced by:

- Heel Stiffness
- Anterior-Posterior positioning of device
- Dorsi-Plantarflexion
- Shoe Characteristics

Consider the following actions if needed:

Symptoms

- Device comes to flat position too early (patient feels he/she is sinking into a hole)
- Climbing over the toe requires extra energy
- Toe feels too stiff
- Knee hyperextends

Action

- Check Wedge instructions (if applicable)
- Shift socket anterior (or device posterior)
- Consider dorsiflexion
- Check heel of the shoe and shoe performance

Symptoms

- Rapid heel to toe movement
- Poor control over prosthesis at initial contact
- Minimal energy return feeling
- Too little push off from the toe
- Knee becomes unstable

Action

- Check Wedge instructions (if applicable)
- Shift socket posterior (or device anterior)
- Consider plantarflexion
- Check heel of the shoe and shoe performance

Heel Wedges (Fig. 8)

The small, medium and large Wedges are used to change heel stiffness. The Wedges can be trimmed using sharp scissors to customize stiffness.

For temporary Wedge placement, use tape to secure the wedge in position.

For permanent Wedge placement

- Roughen the upper and lower surface of the heel with abrasive paper.
- Apply instant adhesive on the lower side of the Wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- If the foot has a split toe, remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon blades.

For removal, the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Noise from Foot

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to two million load cycles.

Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - Label		
Category	Weight (Kg) Moderate Impact	Label Text
1	52	ISO 10328-P3*-52 kg
2	59	ISO 10328-P3*-59 kg
3	68	ISO 10328-P3-68 kg
4	77	ISO 10328-P3-77 kg
5	88	ISO 10328-P4-88 kg
6	100	ISO 10328-P5-100 kg
7	116	ISO 10328-P6-116 kg
8	130	ISO 10328-P7-130 kg
9	147	ISO 10328-P8-147 kg
*Actual test loads reflect maximum body mass		

ISO 10328 - "P" - "m"kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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