
ENGLISH

MD Medical Device

INTENDED USE

The device is intended as a part of a system that replaces function of a missing limb. It connects the knee joint with the proximal part of the prosthesis and allows for inward and outward rotation of the knee joint. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional. The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb amputation and/or congenital deficiency
- No known contraindications

The weight limit for the device is 125 kg.

The device is for moderate impact use, e.g., walking.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Össur replacement parts for the device for ensured durability. The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: For secure fastening use medium strength threadlocker and torque twice to 15 Nm.

Caution: The device must be mounted with the release button in a medial or lateral position and installed directly above the knee joint.

Caution: When using the device in combination with the Unity® vacuum suspension, ensure sufficient tube length, and instruct the user accordingly to avoid damage due to over-stretching or twisting.

USAGE

Warning: Load the prosthesis in the locked state only.

- Press the release button (A) while the prosthesis is not loaded to enable rotation (**Fig. 1**).
- Rotate the device as needed. It can be rotated by 360°.
- The device locks automatically when rotated back to the starting position.

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed. No contact with salt water or chlorinated water is allowed. Dry with a cloth after contact with fresh water or humidity. Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY


Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

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