



Instructions for Use

LOCKING KNEE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	16
IT	Istruzioni per l'uso	20
NO	Bruksanvisning	24
DA	Brugsanvisning	27
SV	Bruksanvisning	31
EL	Οδηγίες Χρήσης	35
FI	Käyttöohjeet	39
NL	Gebruiksaanwijzing	43
PT	Instruções de Utilização	47
PL	Instrukcja użytkowania	51
TR	Kullanım Talimatları	55
RU	Инструкция по использованию	59
JA	取扱説明書	63
ZH	中文说明书	67
KO	사용 설명서	70

DESCRIPTION

The device is a single axis knee joint with a locking mechanism that can be unlocked by means of a Nylon Cord (**Fig. 4: E**). The device remains locked while in a fully extended position but can be unlocked when the patient needs to flex the knee.

The device has a Distal Tube Connector and a Proximal Integrated Male Pyramid Adapter.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low impact use, e.g., gentle walking.

The weight limit for the device is 125 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions. The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Mounting the Manual Release (Fig. 4)

The Nylon Cord (E) should run along the lateral or anterior side of the prosthesis.

Anterior side assembly:

Warning: Failure to adhere to the following instructions may cause the nylon cord to intervene with the locking function of the knee, which may lead to a fall.

1. Mould the Screw Insert (A) into the socket at a desired location.
2. Attach the device to the socket.
3. Fasten the Handle Assembly (B) to the socket with the Handle Fixation Screw (C) and Collar Washer (D).
4. Thread the Nylon Cord (E) through the Cable Guide (F) and Metal Fastener (G).

5. Thread the Nylon Cord (E) through the central hole in the device and again through the Metal Fastener (G).
6. Shorten the Nylon Cord (E) to the length that suits the patient and make sure that the Nylon Cord is slightly tensed.
7. Crimp the Metal Fastener (G) by pliers or pinch the center with a punch and mallet.
8. Secure the Cable Guide (F) to the socket with a small amount of super glue or other means of preference.

Lateral side assembly:

Follow the same instructions as for anterior side assembly except:

- when running the Nylon Cord (E) through the central hole in the device, use the Crimp (H).

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 1)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through midtube (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through midtube (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i.e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 2**).

Tighten the screws with the following torque:

- Tube Clamp Screw (**Fig. 3: A**): 10 Nm

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.
Ensure the patient is familiar with the manual lock release function.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Lubricate the device at least every 6 months with general purpose oil.

1. Flex the knee joint.
2. Place few drops of oil on the joint.
3. Move the knee several times and wipe off excessive oil with a soft cloth.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P6 - 125 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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