







Instructions for Use

OHP3 KNEE

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	16
ES	Instrucciones para el uso	22
IT	Istruzioni per l'uso	28
NO	Bruksanvisning	34
DA	Brugsanvisning	39
SV	Bruksanvisning	44
EL	Οδηγίες Χρήσης	49
FI	Käyttöohjeet	55
NL	Gebruiksaanwijzing	60
PT	Instruções de Utilização	66
PL	Instrukcja użytkowania	71
TR	Kullanım Talimatları	77
RU	Инструкция по использованию	82
JA	取扱説明書	88
ZH	中文说明书	93
KO	사용 설명서	98

DESCRIPTION

The device is a pneumatic polycentric knee joint with a geometric locking system and manual swing phase adjustment.

The device has a distal tube connector and is supplied with one of these four proximal adapter options (**Fig. 1**). Only replace the adapter with one of these adapter options:

- Pyramid Adapter (A)
- 36 mm Thread (B) for use with 3-Prong Adapter
- Loop Adapter (C)
- IKF Adapter (D)

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate impact use, e.g., walking.

The weight limit for the device is 125 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions. The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

36 mm Thread for use with 3-Prong Socket Adapter

Please refer to 3-Prong Socket Adapter Instructions for use.

After lamination, screw the threaded adapter (**Fig. 1: B**) completely into the 3-prong adapter and tighten the lock screw to 10 Nm. After adjustments, the lock screw must be secured with medium strength threadlocker.

Caution: When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without

damaging the adapter so that adapter threads can be fully engaged when mounted on the device.

Caution: When carrying out any later adjustments, do not unscrew the threaded adapter more than one full turn. Threaded adapter might break.

Loop Adapter (Fig. 1)

Lamination of adapter to socket:

- Create a dummy disc to replace the adapter disc during lamination. Fix it to the loop adapter (E).
- Laminate the loop adapter (E) to the socket.
- Protect adapter threads during lamination.
- Remove the dummy disc and attach adapter disc (F) before attaching the loop adapter to the knee.

Caution: Fill the central thread of the loop adapter completely with Plastilin, so that no resin can run into the thread during lamination and the central screw can later be completely screwed in.

Caution: Check that the fine serrations on the inner surfaces of the adjustment discs face towards each other.

All Adapters (Fig. 1)

1. Fit the Washers (L and M) onto the Central Screw (I, J or K),
 - match the conical surfaces together.
 - match the fine serrations surfaces together.
2. Mount the adapter on the knee with the Central Screw, using medium strength threadlocker.
3. Tighten the screw to 35 Nm.

Warning: Screws for other adapters may have different lengths. Do not mix the central screws. Check screw for full thread engagement.

Caution: The IKF adapter has reference marks A and P for Anterior and Posterior (G). Make sure to position correctly.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 2)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through midtube (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through midtube (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).

6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 3**).

Tighten the screws with the following torque:

- Central Screw (**Fig. 1: I, J, or K**): 35 Nm
- Tube Clamp Screw (**Fig. 4: A**): 16 Nm

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

Ensure the patient can activate and deactivate the geometric lock.

Adjustment of stance phase stability (**Fig. 5 and 6**)

Stance phase stability can be adjusted by shifting the knee or by using the supplied wedges to change axle geometry according to the following:

To increase stance phase stability:

- shift knee backwards compared to socket (5A) or
- remove wedge (6B)

To allow for a more dynamic gait:

- shift knee forward compared to socket (5B) or
- Add wedge A (6B) or
- Add wedge B (most dynamic)

To change or remove the wedge (6B), fully flex the knee and loosen the screws (6A). After changing the wedge, apply medium strength threadlocker to the screws and tighten by hand.

For IKF adapter only: Stance flexion adjustment

The back bumper (**Fig. 1: H**) can be changed to increase or decrease stance flexion:

Back Bumper (with holes)		
Yellow	Soft	More stance phase damping
Red	Medium	Standard
Brown	Hard	Less stance phase damping

To replace a bumper, unscrew the central screw. Replace the bumper (holes of back bumper facing up) and reassemble according to assembly instructions above (torque to 35 Nm).

Swing Control (Fig. 7)

The pneumatic system is supplied in its basic setting.

In the basic setting, the Flexion and Extension Valves are opened by 2 full turns. The device can be adjusted to the individual gait.

The Flexion and Extension Valves are labeled in English on the device:

- "Flexion": marks the Flexion valve,
- "Extension": marks the Extension valve.

Adjusting the Flexion and Extension Valves:

The following valves can be adjusted:

- The Flexion Valve: Affects swing flexion resistance.
- The Extension Valve: Affects swing extension resistance.

Start with adjusting the flexion with the patient walking slowly and at faster walking speeds.

For each adjustment, turn the valve in small increments of about a 1/4 of a turn. Always check the result right after each adjustment step.

Possible observations and actions to take are outlined in the following:

- If the patient is walking fast and you observe excessive heel rise, increase the swing flexion resistance by turning the Flexion Valve to the right until the heel rise is normalized.

Afterwards, adjust the extension to harmonize the gait.

- Increase the swing extension resistance by turning the Extension Valve to the right until terminal impact at full extension is reduced. Patient should feel a slight bump at full extension.
- Fine tune the Flexion and Extension Valves until a smooth and secure gait pattern is achieved for slow and fast walking speeds.

Caution: The valves should never be fully closed. Over tightening valves or flexing knee with all valves completely closed may damage valves.

Warning: Valves adjustment must always allow for complete extension of the knee. Excessive swing extension resistance will prevent knee from moving into full extension and may compromise stance security.

USAGE

Knee cap removal and mounting

To remove the knee cap, pull the knee cap straight to the front (**Fig. 8**). When mounting the knee cap, push the knee cap straight on from the front.

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -10°C to 40°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Lubricate piston rod with silicone spray or similar.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P6 - 125 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use
see manufacturer's written instructions on
intended use!

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