




Instructions for Use

OP4 KNEE

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	28
DA	Brugsanvisning	32
SV	Bruksanvisning	36
EL	Οδηγίες Χρήσης	40
FI	Käyttöohjeet	45
NL	Gebruiksaanwijzing	49
PT	Instruções de Utilização	55
PL	Instrukcja użytkowania	60
TR	Kullanım Talimatları	65
RU	Инструкция по использованию	69
JA	取扱説明書	74
ZH	中文说明书	78
KO	사용 설명서	82

DESCRIPTION

The device is a pneumatic single axis knee joint with a weight-activated friction brake, and an individually adjustable swing phase.

The device has a Distal Tube Connector and a Proximal Integrated Male Pyramid Adapter.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to moderate impact use, e.g., walking.

The weight limit for the device is 100 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions. The device is for single patient use.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 1)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through midtube (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through midtube (A)

4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 2**).

Tighten the screws with the following torque:

- Tube Clamp Screw (**Fig. 3: A**): 16 Nm

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

- Ensure the patient is familiar with the functioning of the device.
- Ensure the patient is familiar with the weight-activated brake function.

Stance phase stability (Fig. 4):

Start carefully with the device's brake function with the patient at the dynamic alignment and check whether this basic setting matches the patient's needs. If not, please proceed as follows:

Adjustment using the Load Regulating Screw (A)

The Load Regulating Screw can be used to adjust:

- the point at which the brake acts when (weight-) loading starts,
- and the response of the brake during the dynamic alignment.

Variant 1 - Brake effect shall be induced later:

Turn the Load Regulating Screw (A) to the right which increases the weight needed by the patient to initiate the brake function.

Variant 2 - Brake effect shall be induced earlier:

Turn the Load Regulating Screw (A) to the left which reduces the weight needed by the patient to initiate the brake function.

Note: Always adjust the Load Regulating Screw (A) carefully. Note that you can adjust it against the slight resistance of the threadlocker.

Note: The Load Regulating Screw (A) must not be screwed beyond the device's brake body. Otherwise, the function of the device cannot be guaranteed.

Readjusting the brake unit after prolonged use

After a prolonged period of use, play may develop in the device which can be readjusted using the Service Screw (B).

Turn the Service Screw by small steps to the right.

Carefully check whether the brake function of the device now suits the patient's individual needs. If not, readjust the Service Screw again until the correct setting has been found.

Note: Always adjust the Service Screw (B) carefully. Note that you can adjust it against the slight resistance of the thread locker.

Caution: flexion and extension must be possible at all settings.

Restoring the basic setting

If necessary, the basic setting of the device can be reset as follows:

Step 1: Carefully turn the Load Regulating Screw (4A) to the right until a slight resistance is felt.

Step 2: Then turn the Load Regulating Screw three turns back to the left to restore the basic setting.

Swing Control (Fig. 5)

The pneumatic system is supplied in its basic setting.

In the basic setting, the Flexion and Extension Valves are opened by 2 full turns. The device can be adjusted to the individual gait.

The Flexion and Extension Valves are labeled in English on the device:

- "Flexion": marks the Flexion valve,
- "Extension": marks the Extension valve.

Adjusting the Flexion and Extension Valves:

The following valves can be adjusted:

- The Flexion Valve: Affects swing flexion resistance.
- The Extension Valve: Affects swing extension resistance.

Start with adjusting the flexion with the patient walking slowly and at faster walking speeds.

For each adjustment, turn the valve in small increments of about a ¼ of a turn. Always check the result right after each adjustment step.

Possible observations and actions to take are outlined in the following:

- If the patient is walking fast and you observe excessive heel rise, increase the swing flexion resistance by turning the Flexion Valve to the right until the heel rise is normalized.

Afterwards, adjust the extension to harmonize the gait.

- Increase the swing extension resistance by turning the Extension Valve to the right until terminal impact at full extension is reduced. Patient should feel a slight bump at full extension.
- Fine tune the Flexion and Extension Valves until a smooth and secure gait pattern is achieved for slow and fast walking speeds.

Caution: The valves should never be fully closed. Over tightening valves or flexing knee with all valves completely closed may damage valves.

Note: Flexion and extension must be possible at all settings.

Note: It may be necessary to readjust the brake again after setting the swing phase.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

Check for damage, excessive wear and dirt.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Lubricate piston rod with silicone spray or similar.

Caution: Do not use talcum powder for lubrication.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P5 - 100 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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