Rev. 3 - 05/12_103649

**Bauerfeind AG**

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EN Dear client,

Please read these instructions for use carefully. If you have any questions, please contact our customer service department. MOS Genu® is a therapeutic device designed to aid your treatment. It helps to protect the knee during therapy and assists muscle development during physical therapy exercises.

Indications

Non-surgical

- Cruciate and lateral ligament injuries
- Complex instabilities including genu recurvatum
- Support for the affected compartment in bowleg/knock-knee arthritis

Post-op. (short MOS Genu)

- Operations on the ligaments, i.e. artificial ligaments
- Meniscal suturing and meniscal replacement (meniscal implantation)

Post-op. (long MOS Genu)

- After axial correction of the shinbone (tibial osteotomy)
- After complex reconstructions of the ligaments
- In individual cases it is also indicated for fractures located close to joints

Side effects

Side effects involving the body as a whole have not been reported side date. Correct fitting is assumed. Any supports and orthoses¹ applied externally to the body can, if tightened excessively, lead to local pressure signs or, rarely, constrict the underlying blood vessels or nerves.

Contraindications

Hypersensitive reactions harmful to health have not been reported to date. In the following conditions such aids should only be fitted and worn after consultation with your doctor:

1. Skin disorders/injuries in the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness and excessive heat build-up.
2. Varicose veins.
3. Impaired sensation and circulatory disorders of the legs, e.g. diabetes mellitus.
4. Impaired lymph drainage – including soft tissue swellings of uncertain origin located remotely from the fitted aid.

Determining the Size

Please determine the correct size according to the details given on the back of the packaging.

Procedure for First Time Application

1. Preparation

- Loosen the 4 adjustment cams above and below the knee ring.
- Mark the theoretical anatomical knee joint centre at 30° flexion.
- Please choose the size of the knee joint padding for the individual. MOS Genu comes ready-fitted with post-op pads on the inner side of the joint for post-operative use.

2. Static fitting

IMPORTANT!

The shell sections and the straps must lie flat on the skin. Do not use sharp-edged tools to contour the orthosis. Most adjustments can be performed manually.

- The initial fitting is carried out without the knee pads.
- The MOS Genu is placed against the knee joint in 30 degree flexion such that the mechanical pivot point of the MOS Genu corresponds with Nietert's compromise axis.
- Check that the support shells lie flat against the knee both vertically and circumferentially; adjust if necessary.
- A round adjusting hook with safety cap should be used for major adjustments.
- Select the knee pads of the required thickness, fasten and turn the Velcro-backed pads until flat contact with the knee is achieved. The post-op pads must be removed before the knee padding is attached using the Velcro straps.
- Centre the distal (lower) popliteal fossa strap with the integrated knee pad, tighten and secure on both sides.
- Take the proximal (upper) popliteal fossa strap and proceed as for the distal strap.
- Now attach, centre and secure the two broader distal and proximal (bottom and top) straps and, on the long version, the straps of the extension shells.
- The brace is held in the individual bowleg/knock-knee position whilst the straps are attached. The adjusting cams should now be screwed tight to maintain this position. Important! Do not attempt to move adjusting cams once tightened, otherwise the material will be stressed.
- On the long version, the height of the distal (lower) and proximal (upper) extension shells can be adjusted to the full extent of the lever arm.
- If necessary, the protective shinbone padding, also supplied, can be affixed in order to relieve pressure pain in the shinbone area - or if swelling is present in the area of the tibial tubercle.

3. Adjustment examples

- Two stop pins can be set to limit joint flexion and extension.
- The extent of permitted movement can be read off, in degrees, in the window on the outer aspect of the knee.
- Flexion/extension can be restricted by placing the stop pins at the appropriate positions.
- Angles of 0°, 20°, 40°, 60° etc. can be set with the stop pin in the long slot and angles of 10°, 30°, 50° etc. with the stop pin in the short slot.
- The orthosis is shipped with a 0°/0°/90° flexion/extension initial setting.


Fig. 1: Schematic diagram of possible joint settings

4. Dynamic fitting

- Check the effect of MOS Genu on the patient's movement.
- For gait anomalies, loosen the adjustment cams and in 20 steps let the orthosis run in, then retighten adjustment cams. Important! Do not attempt to move adjusting cams once tightened, otherwise the material will be stressed.
- Repeat this procedure until the desired gait is reached.
- TORX screws have been used exclusively on the MOS Genu. A matching screwdriver is also supplied. TORX screws allow for high tightening torques as it is virtually impossible for the screwdriver to slip out of the TORX drive.

5. Bowleg/knock-knee settings

- The graduated adjusting cams allow the brace to be placed in the precise bowleg/knock-knee position of the leg.
- Loosen the cam screws and set the guide cams first (identified by spots in the illustration). The follower cams will turn in the opposite direction.
- Precise adjustment is achieved by alternating the rotation of both cams distally (lower) and proximally (upper).

¹ Orthosis = orthopaedic appliance used to stabilize, relieve, immobilize, control or correct the limbs or the torso

- The distal (bottom) and proximal (top) scale settings on the cam pairs should be the same. Adding the distal (bottom) and proximal (top) readings together will give the individual setting.
- In appropriate applications, targeted relief of the ligaments on the inner or outer aspects of the knee can be achieved by correcting the bowleg/knock-knee setting.




deformity	upper (proximal)		lower (distal)		setting
					
Genu-Varus		3°	+	4°	$\Delta = 7^\circ$

Fig. 2: Varus setting of 7° overall on the left leg

6. Dismantling the long version

- Depending on the application and the progress made by the patient in the rehabilitation phase, it is possible to simply dismantle the brace and turn it into the short version.
- Release and remove the screws on the distal (lower) and proximal (upper) support shells.
- Pull the distal (lower) and proximal (upper) extension shells out of the support shells.
- Compensate for the resulting hollows on the inner aspects of the shells with the padded strips provided.
- Place the orthosis in position and check that the fit is correct.

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- A 4-point system is implemented with the aid of the two popliteal fossa pads (Fig. 3), e.g. for an anterior instability (drawer sign).
- For a posterior instability (drawer sign) or back-knee (genu recurvatum), there is a long popliteal cross-strap available to replace the short popliteal straps.
- Twist the cross-strap through 180° and thread through the pad loops. Pull through the fixing loops on the distal (lower) joint ring and cross over in a proximal (upper) direction and secure (Fig. 4).
- The popliteal cross-strap can also be used after surgery to make the brace more comfortable to wear.


Fig. 3: ACL Treatment



Fig. 4: PCL Treatment

- Changes must not be made for the product. Non compliance can reduce the performance of the product and thus no liability will be accepted.
- If you notice any unusual changes in the patient (e.g. increased symptoms), he or she must consult their doctor without delay.
- A combination with other products (e.g. compression stockings) must be discussed beforehand with your doctor.
- Avoid the use of ointments, lotions or any substances including grease or acids.
- Spare parts and accessories are available.
- MOS Genu is not designed to prevent injuries, especially when engaging contact sports. Contact sports inherently pose a greater risk of injury.
- The product can be disposed of without any restrictions.

Product care

⚠ Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!

The aluminium parts of the support are coated and can be cleaned with a pH-neutral soap. We recommend washing the pads at 30° C on a gentle wash cycle and hand-washing the straps at 30° C.

This product has been examined by our own quality control system. However, should you have any complaint, please contact your retailer or our customer service department.

Version: January 2011

⚠ Important information

- MOS Genu is a prescribed product that should only be worn under a physician's guidance. MOS Genu must be worn in accordance with the information contained in these instructions for use and for the indications listed only.
- The product is intended to be fitted only once and to one patient.
- No liability is accepted for improper use of the product.

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Bauerfeind AG die Konformität zur